



PMEL Antibody, FITC conjugated

Product Code	CSB-PA021324LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40967
Immunogen	Recombinant Human Melanocyte protein PMEL protein (25-467AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Melanocyte protein PMEL (ME20-M) (ME20M) (Melanocyte protein Pmel 17) (Melanocytes lineage-specific antigen GP100) (Melanoma-associated ME20 antigen) (P1) (P100) (Premelanosome protein) (Silver locus protein homolog) [Cleaved into: M-alpha (95 kDa melanocyte-specific secreted glycoprotein) (P26) (Secreted melanoma-associated ME20 antigen) (ME20-S) (ME20S); M-beta], PMEL, D12S53E PMEL17 SILV
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	PMEL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.