



gH Antibody, Biotin conjugated

Product Code	CSB-PA319015LD01HWV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12824
Immunogen	Recombinant Human cytomegalovirus Envelope glycoprotein H protein (24-195AA)
Raised In	Rabbit
Species Reactivity	Human cytomegalovirus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Envelope glycoprotein H (gH), gH, UL75
Immunogen Species	Human cytomegalovirus
Research Area	Others
Target Names	gH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.