



NPR1 Antibody, HRP conjugated

Product Code	CSB-PA016023LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16066
Immunogen	Recombinant Human Atrial natriuretic peptide receptor 1 protein (33-473AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Atrial natriuretic peptide receptor 1 (EC 4.6.1.2) (Atrial natriuretic peptide receptor type A) (ANP-A) (ANPR-A) (NPR-A) (Guanylate cyclase A) (GC-A), NPR1, ANPRA
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	NPR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.