



thiO Antibody, FITC conjugated

Product Code	CSB-PA521738LC01BRJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O31616
Immunogen	Recombinant Bacillus subtilis Glycine oxidase protein (1-369AA)
Raised In	Rabbit
Species Reactivity	Bacillus subtilis
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glycine oxidase (GO) (EC 1.4.3.19), thiO, goxB yjbR
Immunogen Species	Bacillus subtilis (strain 168)
Research Area	Others
Target Names	thiO
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.