



ALB Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA001561LC01MOV |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | A2V9Z4 |
| Immunogen | Recombinant Macaca fascicularis Serum albumin protein (25-608AA) |
| Raised In | Rabbit |
| Species Reactivity | Macaca fascicularis |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Serum albumin, ALB |
| Immunogen Species | Macaca fascicularis (Crab-eating macaque) (Cynomolgus monkey) |
| Research Area | Others |
| Target Names | ALB |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |