



ALB Antibody, Biotin conjugated

Product Code	CSB-PA001561LD01MOV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A2V9Z4
Immunogen	Recombinant Macaca fascicularis Serum albumin protein (25-608AA)
Raised In	Rabbit
Species Reactivity	Macaca fascicularis
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serum albumin, ALB
Immunogen Species	Macaca fascicularis (Crab-eating macaque) (Cynomolgus monkey)
Research Area	Others
Target Names	ALB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.