



OV16 Antibody, Biotin conjugated

Product Code	CSB-PA341251HD01OEG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31729
Immunogen	Recombinant Onchocerca volvulus OV-16 antigen protein (17-197AA)
Raised In	Rabbit
Species Reactivity	Onchocerca volvulus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	OV-16 antigen, OV16
Immunogen Species	Onchocerca volvulus
Research Area	Others
Target Names	OV16
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.