




Efnb3 Antibody

Product Code	CSB-PA007467LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O35393
Immunogen	Recombinant Mouse Ephrin-B3 protein (28-227AA)
Raised In	Rabbit
Species Reactivity	Mouse, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ephrin-B3, Efnb3
Immunogen Species	Mus musculus (Mouse)
Research Area	Cardiovascular
Target Names	Efnb3

Image



Western blot
 All lanes: Efnb3 antibody at 12µg/ml + HepG2 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 36 kDa
 Observed band size: 36 kDa

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.