



Adipoq Antibody, HRP conjugated

Product Code	CSB-PA07946B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q60994
Immunogen	Recombinant Mouse Adiponectin protein (18-247AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adiponectin (30 kDa adipocyte complement-related protein) (Adipocyte complement-related 30 kDa protein) (ACRP30) (Adipocyte, C1q and collagen domain-containing protein) (Adipocyte-specific protein AdipoQ), Adipoq, Acdc Acrp30 Apm1
Immunogen Species	Mus musculus (Mouse)
Research Area	Cardiovascular
Target Names	Adipoq
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.