



# Adipoq Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA07946B0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q60994
<b>Immunogen</b>	Recombinant Mouse Adiponectin protein (18-247AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Adiponectin (30 kDa adipocyte complement-related protein) (Adipocyte complement-related 30 kDa protein) (ACRP30) (Adipocyte, C1q and collagen domain-containing protein) (Adipocyte-specific protein AdipoQ), Adipoq, Acdc Acrp30 Apm1
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Research Area</b>	Cardiovascular
<b>Target Names</b>	Adipoq