



HIF1A Antibody, HRP conjugated

Product Code	CSB-PA12807B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16665
Immunogen	Recombinant Human Hypoxia-inducible factor 1-alpha protein (579-826AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hypoxia-inducible factor 1-alpha (HIF-1-alpha) (HIF1-alpha) (ARNT-interacting protein) (Basic-helix-loop-helix-PAS protein MOP1) (Class E basic helix-loop-helix protein 78) (bHLHe78) (Member of PAS protein 1) (PAS domain-containing protein 8), HIF1A, BHLHE78 MOP1 PASD8
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	HIF1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.