



Kif1c Antibody

Product Code	CSB-PA012320LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O35787
Immunogen	Recombinant Rat Kinesin-like protein KIF1C protein (1-272AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kinesin-like protein KIF1C (Kinesin-like protein KIF1D), Kif1c, Kif1d
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Immunology
Target Names	Kif1c