



trxA Antibody, Biotin conjugated

Product Code	CSB-PA18667D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0AA25
Immunogen	Recombinant Escherichia coli Thioredoxin 1 protein (2-109AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thioredoxin 1 (Trx-1), trxA, fipA tsnC
Immunogen Species	Escherichia coli
Research Area	Signal Transduction
Target Names	trxA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.