



Tpsb2 Antibody

Product Code	CSB-PA024128LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50343
Immunogen	Recombinant Rat Tryptase beta-2 protein (30-274AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tryptase beta-2 (Tryptase-2) (EC 3.4.21.59) (Mast cell protease 6) (rMCP-6), Tpsb2, Mcp6 Mcpt6
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Immunology
Target Names	Tpsb2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.