



TRBC1 Antibody, Biotin conjugated

Product Code	CSB-PA024291LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01850
Immunogen	Recombinant Human T cell receptor beta constant 1 protein (1-176AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-cell receptor beta-1 chain C region, TRBC1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TRBC1