



ssb Antibody, HRP conjugated

Product Code	CSB-PA770616LB01BQE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q81JI3
Immunogen	Recombinant Bacillus anthracis Single-stranded DNA-binding protein (1-172AA)
Raised In	Rabbit
Species Reactivity	Bacillus anthracis
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Single-stranded DNA-binding protein (SSB), ssb, ssb-1
Immunogen Species	Bacillus anthracis
Research Area	Others
Target Names	ssb
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.