



GALE Antibody, Biotin conjugated

Product Code	CSB-PA613499HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14376
Immunogen	Recombinant Human UDP-glucose 4-epimerase protein (1-348AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UDP-glucose 4-epimerase (EC 5.1.3.2) (Galactowaldenase) (UDP-N-acetylgalactosamine 4-epimerase) (UDP-GalNAc 4-epimerase) (UDP-N-acetylglucosamine 4-epimerase) (UDP-GlcNAc 4-epimerase) (EC 5.1.3.7) (UDP-galactose 4-epimerase), GALE
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GALE
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.