



Genome polyprotein Antibody, HRP conjugated

Product Code	CSB-PA333180LB01HFD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26664
Immunogen	Recombinant Hepatitis C virus genotype 1a Genome polyprotein protein (192-325AA)
Raised In	Rabbit
Species Reactivity	Hepatitis C virus genotype 1a
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Genome polyprotein [Cleaved into: Core protein p21 (Capsid protein C) (p21); Core protein p19; Envelope glycoprotein E1 (gp32) (gp35); Envelope glycoprotein E2 (NS1) (gp68) (gp70); p7; Protease NS2-3 (p23) (EC 3.4.22.-); Serine protease NS3 (EC 3.4.21.98) (EC 3.6.1.15) (EC 3.6.4.13) (Hepacivirin) (NS3P) (p70); Non-structural protein 4A (NS4A) (p8); Non-structural protein 4B (NS4B) (p27); Non-structural protein 5A (NS5A) (p56); RNA-directed RNA polymerase (EC 2.7.7.48) (NS5B) (p68)]
Immunogen Species	Hepatitis C virus genotype 1a
Research Area	Microbiology