



# Genome polyprotein Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA333180LD01HFD
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P26664
<b>Immunogen</b>	Recombinant Hepatitis C virus genotype 1a Genome polyprotein protein (192-325AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Hepatitis C virus genotype 1a
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Genome polyprotein [Cleaved into: Core protein p21 (Capsid protein C) (p21); Core protein p19; Envelope glycoprotein E1 (gp32) (gp35); Envelope glycoprotein E2 (NS1) (gp68) (gp70); p7; Protease NS2-3 (p23) (EC 3.4.22.-); Serine protease NS3 (EC 3.4.21.98) (EC 3.6.1.15) (EC 3.6.4.13) (Hepacivirin) (NS3P) (p70); Non-structural protein 4A (NS4A) (p8); Non-structural protein 4B (NS4B) (p27); Non-structural protein 5A (NS5A) (p56); RNA-directed RNA polymerase (EC 2.7.7.48) (NS5B) (p68)]
<b>Immunogen Species</b>	Hepatitis C virus genotype 1a
<b>Research Area</b>	Microbiology
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.