



PRR3 Antibody

Product Code	CSB-PA018794LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P79522
Immunogen	Recombinant Human Proline-rich protein 3 protein (1-188AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proline-rich protein 3 (MHC class I region proline-rich protein CAT56), PRR3, CAT56
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PRR3