



CEP83 Antibody

Product Code	CSB-PA896730LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y592
Immunogen	Recombinant Human Centrosomal protein of 83 kDa protein (1-568AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Centrosomal protein of 83 kDa (Cep83) (Coiled-coil domain-containing protein 41) (Renal carcinoma antigen NY-REN-58), CEP83, CCDC41
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	CEP83
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.