




ALOX15B Antibody

Product Code	CSB-PA001622LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15296
Immunogen	Recombinant Human Arachidonate 15-lipoxygenase B protein (1-676AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IF; Recommended dilution:IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Arachidonate 15-lipoxygenase B (15-LOX-B) (EC 1.13.11.33) (15-lipoxygenase 2) (15-LOX-2) (Arachidonate 15-lipoxygenase type II) (Linoleate 13-lipoxygenase 15-LOb) (EC 1.13.11.-), ALOX15B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ALOX15B

Image



Immunofluorescence staining of A549 cells with CSB-PA001622LA01HU at 1:166, counter-stained with DAPI. The cells were fixed in 4% formaldehyde, permeabilized using 0.2% Triton X-100 and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.