




# PLEKHB2 Antibody

<b>Product Code</b>	CSB-PA846599LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q96CS7
<b>Immunogen</b>	Recombinant Human Pleckstrin homology domain-containing family B member 2 protein (1-222AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Pleckstrin homology domain-containing family B member 2 (PH domain-containing family B member 2) (Evectin-2), PLEKHB2, EVT2
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	PLEKHB2

## Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA846599LA01HU at dilution of 1:100