



IFNA1 Antibody, HRP conjugated

Product Code	CSB-PA365592HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01562
Immunogen	Recombinant Human Interferon alpha-1/13 protein (24-189AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon alpha-1/13 (IFN-alpha-1/13) (Interferon alpha-D) (LeIF D), IFNA1; IFNA13
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IFNA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.