



WDR5 Antibody

Product Code	CSB-PA026038HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61964
Immunogen	Recombinant Human WD repeat-containing protein 5 protein (2-334AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	WD repeat-containing protein 5 (BMP2-induced 3-kb gene protein), WDR5, BIG3
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	WDR5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.