



NTS Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA016136HC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P30990 |
| Immunogen | Recombinant Human Neurotensin/neuromedin N protein (24-143AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Neurotensin/neuromedin N [Cleaved into: Large neuromedin N (NmN-125); Neuromedin N (NN) (NmN); Neurotensin (NT); Tail peptide], NTS |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Neuroscience |
| Target Names | NTS |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |