



gluD Antibody

Product Code	CSB-PA335284HA01DUN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27346
Immunogen	Recombinant Clostridium difficile NAD-specific glutamate dehydrogenase protein (1-421AA)
Raised In	Rabbit
Species Reactivity	Clostridium difficile
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NAD-specific glutamate dehydrogenase (NAD-GDH) (EC 1.4.1.2), gluD
Immunogen Species	Clostridioides difficile (Peptoclostridium difficile)
Research Area	Others
Target Names	gluD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.