



# Recombinant Mouse Pyruvate dehydrogenase E1 component subunit alpha, testis-specific form, mitochondrial (Pdha2)

<b>Product Code</b>	CSB-BP017716MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P35487
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	FSNDATCDIK KCDLYRLEEG PPTSTVLTRA EALKYYRTMQ VIRRMELKAD QLYKQKFIRG FCHLCDGQEA CCVGLEAGIN PTDHVITSYR AHGFCYTRGL SVKSILAELELT GRKGGCAKGGK GGSMHMYGKN FYGGNGIVGA QVPLGAGVAF ACKYLKNGQV CLALYGDGAA NQGQVFEAYN MSALWKLPCV FICENNL YGM GTSNERSAAS TDYHKKGFII PGLRVNGMDI LCVREATKFA ADHCRSGKGP IVMELQTYRY HGHSMSDPGI SYRSREEVHN VRSKSDPIML LRERIISNNL SNIEELKEID ADVKKEVEDA AQFATTDPEP AVEDIANYLY HQDPPFEVRG AHKWLKYKSH S
<b>Source</b>	Baculovirus
<b>Target Names</b>	Pdha2
<b>Protein Names</b>	Recommended name: Pyruvate dehydrogenase E1 component subunit alpha, testis-specific form, mitochondrial EC= 1.2.4.1 Alternative name(s): PDHE1-A type II
<b>Expression Region</b>	31-391
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.