



scalM Antibody, FITC conjugated

Product Code	CSB-PA528741LC01FNZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O52692
Immunogen	Recombinant Streptomyces caespitosus Modification methylase Scal protein (1-304AA)
Raised In	Rabbit
Species Reactivity	Streptomyces caespitosus
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Modification methylase Scal (M.Scal) (EC 2.1.1.113) (N-4 cytosine-specific methyltransferase Scal), scalM
Immunogen Species	Streptomyces caespitosus
Research Area	Others
Target Names	scalM
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.