



GADD45A Antibody

| | |
|---------------------|---|
| Product Code | CSB-PA009161HA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P24522 |
| Immunogen | Recombinant Human Growth arrest and DNA damage-inducible protein GADD45 alpha protein (1-165AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Growth arrest and DNA damage-inducible protein GADD45 alpha (DNA damage-inducible transcript 1 protein) (DDIT-1), GADD45A, DDIT1 GADD45 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling |
| Target Names | GADD45A |