



GADD45A Antibody

Product Code	CSB-PA009161HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24522
Immunogen	Recombinant Human Growth arrest and DNA damage-inducible protein GADD45 alpha protein (1-165AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Growth arrest and DNA damage-inducible protein GADD45 alpha (DNA damage-inducible transcript 1 protein) (DDIT-1), GADD45A, DDIT1 GADD45
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	GADD45A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.