





HAS3 Antibody

Product Code	CSB-PA010142ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00219
Immunogen	Recombinant Human Hyaluronan synthase 3 protein (339-553AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hyaluronan synthase 3 (EC 2.4.1.212) (Hyaluronate synthase 3) (Hyaluronic acid synthase 3) (HA synthase 3), HAS3
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	HAS3


Image



Immunohistochemistry of paraffin-embedded human thyroid tissue using CSB-PA010142ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human cervical cancer using CSB-PA010142ESR1HU at dilution of 1:100



Immunofluorescent analysis of U251 cells using CSB-PA010142ESR1HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.