




NRP1 Antibody

Product Code	CSB-PA016091ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14786
Immunogen	Recombinant Human Neuropilin-1 protein (684-923AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuropilin-1 (Vascular endothelial cell growth factor 165 receptor) (CD antigen CD304), NRP1, NRP VEGF165R
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	NRP1

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA016091ESR1HU at dilution of 1:100