



PAG1 Antibody

Product Code	CSB-PA885746GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NWQ8
Immunogen	Human PAG1
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	Csk-binding protein,Transmembrane adapter protein PAG,Transmembrane phosphoprotein Cbp,CBP, PAG,PAG1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PAG1
Target Details	This protein is a type III transmembrane adaptor protein that binds to the tyrosine kinase csk protein. It is thought to be involved in the regulation of T cell activation.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.