




ALOX15B Antibody

Product Code	CSB-PA001622ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15296
Immunogen	Recombinant Human Arachidonate 15-lipoxygenase B protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Arachidonate 15-lipoxygenase B (15-LOX-B) (EC 1.13.11.33) (15-lipoxygenase 2) (15-LOX-2) (Arachidonate 15-lipoxygenase type II) (Linoleate 13-lipoxygenase 15-LOb) (EC 1.13.11.-), ALOX15B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ALOX15B

Image



Western blot

All lanes: Arachidonate 15-lipoxygenase B antibody at 11µg/ml

Lane 1: MCF-7 whole cell lysate


Lane 2: 293T whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 76, 68, 70, 73 kDa

Observed band size: 76 kDa



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA001622ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.