





PDPK1 Antibody

Product Code	CSB-PA017738ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15530
Immunogen	Recombinant Human 3-phosphoinositide-dependent protein kinase 1 protein (277-556AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	3-phosphoinositide-dependent protein kinase 1 (hPDK1) (EC 2.7.11.1), PDPK1, PDK1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PDPK1


Image



Immunohistochemistry of paraffin-embedded human thyroid tissue using CSB-PA017738ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA017738ESR1HU at dilution of 1:100



Immunofluorescent analysis of MCF-7 cells using CSB-PA017738ESR1HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)