




TRAF3IP2 Antibody

Product Code	CSB-PA024150ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43734
Immunogen	Recombinant Human Adapter protein CIKS protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adapter protein CIKS (Connection to IKK and SAPK/JNK) (Nuclear factor NF-kappa-B activator 1) (ACT1) (TRAF3-interacting protein 2), TRAF3IP2, C6orf2 C6orf4 C6orf5 C6orf6
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TRAF3IP2

Image


Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA024150ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.