





BNIP3L Antibody

Product Code	CSB-PA002767ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60238
Immunogen	Recombinant Human BCL2/adenovirus E1B 19 kDa protein-interacting protein 3-like protein (1-187AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	BCL2/adenovirus E1B 19 kDa protein-interacting protein 3-like (Adenovirus E1B19K-binding protein B5) (BCL2/adenovirus E1B 19 kDa protein-interacting protein 3A) (NIP3-like protein X) (NIP3L), BNIP3L, BNIP3A BNIP3H NIX
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	BNIP3L

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA002767ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA002767ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.