





RNF8 Antibody

Product Code	CSB-PA019898ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O76064
Immunogen	Recombinant Human E3 ubiquitin-protein ligase RNF8 protein (1-260AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase RNF8 (hRNF8) (EC 2.3.2.27) (RING finger protein 8) (RING-type E3 ubiquitin transferase RNF8), RNF8, KIAA0646
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RNF8

Image



Western blot
 All lanes: RNF8 antibody at 11µg/ml
 Lane 1: Hela whole cell lysate
 Lane 2: 293T whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 56, 11, 51 kDa
 Observed band size: 56 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA019898ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.