





BMP15 Antibody

Product Code	CSB-PA002735ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95972
Immunogen	Recombinant Human Bone morphogenetic protein 15 protein (143-392AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Bone morphogenetic protein 15 (BMP-15) (Growth/differentiation factor 9B) (GDF-9B), BMP15, GDF9B
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	BMP15

Image



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA002735ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA002735ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.