





SERPINA3 Antibody

Product Code	CSB-PA021059ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01011
Immunogen	Recombinant Human Alpha-1-antichymotrypsin protein (24-423AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-1-antichymotrypsin (ACT) (Cell growth-inhibiting gene 24/25 protein) (Serpina3) [Cleaved into: Alpha-1-antichymotrypsin His-Pro-less], SERPINA3, AACT
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SERPINA3

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA021059ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human lung tissue using CSB-PA021059ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.