





ARG1 Antibody

Product Code	CSB-PA002005ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05089
Immunogen	Recombinant Human Arginase-1 protein (1-322AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Arginase-1 (EC 3.5.3.1) (Liver-type arginase) (Type I arginase), ARG1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ARG1

Image



Western blot
 All lanes: Arginase-1 antibody at 3µg/ml +
 Mouse liver tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 35, 36, 26 kDa
 Observed band size: 35 kDa



Immunohistochemistry of paraffin-embedded
 human liver tissue using CSB-
 PA002005ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.