





NR2C1 Antibody

Product Code	CSB-PA016051ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13056
Immunogen	Recombinant Human Nuclear receptor subfamily 2 group C member 1 protein (1-240AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear receptor subfamily 2 group C member 1 (Orphan nuclear receptor TR2) (Testicular receptor 2), NR2C1, TR2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NR2C1

Image



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA016051ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human lung tissue using CSB-PA016051ESR2HU at dilution of 1:100