





GRN Antibody

Product Code	CSB-PA009939ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28799
Immunogen	Recombinant Human Granulins protein (20-310AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Granulins (Proepithelin) (PEPI) [Cleaved into: Acrogranin (Glycoprotein of 88 Kda) (GP88) (Glycoprotein 88) (Progranulin); Paragranulin; Granulin-1 (Granulin G); Granulin-2 (Granulin F); Granulin-3 (Granulin B); Granulin-4 (Granulin A); Granulin-5 (Granulin C); Granulin-6 (Granulin D); Granulin-7 (Granulin E)], GRN
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GRN

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA009939ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human melanoma using CSB-PA009939ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.