




GRIA1 Antibody

Product Code	CSB-PA009898ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42261
Immunogen	Recombinant Human Glutamate receptor 1 protein (19-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutamate receptor 1 (GluR-1) (AMPA-selective glutamate receptor 1) (GluR-A) (GluR-K1) (Glutamate receptor ionotropic, AMPA 1) (GluA1), GRIA1, GLUH1, GLUR1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GRIA1

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA009898ESR1HU at dilution of 1:100