




IFNAR2 Antibody

Product Code	CSB-PA011047ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48551
Immunogen	Recombinant Human Interferon alpha/beta receptor 2 protein (27-243AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon alpha/beta receptor 2 (IFN-R-2) (IFN-alpha binding protein) (IFN-alpha/beta receptor 2) (Interferon alpha binding protein) (Type I interferon receptor 2), IFNAR2, IFNABR IFNARB
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IFNAR2

Image



Western blot

All lanes: IFNAR2 antibody at 10µg/ml + MCF-7 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 58, 38, 28 kDa

Observed band size: 58 kDa



Immunohistochemistry of paraffin-embedded human endometrial cancer using CSB-PA011047ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.