




# SULT1A1 Antibody

<b>Product Code</b>	CSB-PA022933ESR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P50225
<b>Immunogen</b>	Recombinant Human Sulfotransferase 1A1 protein (1-295AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Sulfotransferase 1A1 (ST1A1) (EC 2.8.2.1) (Aryl sulfotransferase 1) (HAST1/HAST2) (Phenol sulfotransferase 1) (Phenol-sulfating phenol sulfotransferase 1) (P-PST 1) (ST1A3) (Thermostable phenol sulfotransferase) (Ts-PST), SULT1A1, STP STP1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	SULT1A1

**Image**


Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA022933ESR1HU at dilution of 1:100

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.