




ATXN1 Antibody

Product Code	CSB-PA002438ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54253
Immunogen	Recombinant Human Ataxin-1 protein (586-815AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ataxin-1 (Spinocerebellar ataxia type 1 protein), ATXN1, ATX1 SCA1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ATXN1

Image



Immunohistochemistry of paraffin-embedded human gallbladder tissue using CSB-PA002438ESR2HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.