




COL2A1 Antibody

Product Code	CSB-PA005739ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02458
Immunogen	Recombinant Human Collagen alpha-1(II) chain protein (1228-1487AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Collagen alpha-1(II) chain (Alpha-1 type II collagen) [Cleaved into: Collagen alpha-1(II) chain; Chondrocalcin], COL2A1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	COL2A1

Image



Western blot


All lanes: Collagen alpha-1 (II) chain antibody at 10µg/ml + 293T whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 142, 135, 30 kDa

Observed band size: 142 kDa



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA005739ESR2HU at dilution of 1:100

