




CD58 Antibody

Product Code	CSB-PA004946ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19256
Immunogen	Recombinant Human Lymphocyte function-associated antigen 3 protein (40-220AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lymphocyte function-associated antigen 3 (Ag3) (Surface glycoprotein LFA-3) (CD antigen CD58), CD58, LFA3
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD58

Image



Western blot

All lanes: Lymphocyte function-associated antigen 3 antibody at 2µg/ml

Lane 1: Hela whole cell lysate


Lane 2: Jurkat whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 29, 27, 28 kDa

Observed band size: 29 kDa



Immunohistochemistry of paraffin-embedded human thymus tissue using CSB-PA004946ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.